

Coverlet® Extremely Absorbent Bandages.



Full Range of Shapes & Sizes.

Latex-Free Fabric Dressings.

Now including an Extra Large patch





Extremely Absorbent Bandages.

For maximum absorption and protection, give your patients Coverlet Adhesive Bandages to aid in their healing process.




Coverlet Adhesive Bandages are made to absorb liquid fast; faster than any of the competitors tested.* In addition, every Coverlet features an extra-large wound pad capable of absorbing 10 times its own weight in water. The 360-degree adhesive surrounding the island pad helps seal the afflicted area from dirt and contamination – keeping wounds clean and providing an environment conducive to healing.



Now available in 13 different shapes and sizes, Coverlet is perfect for use on knees, elbows, chins, noses and knuckles. Its elastic fabric moves with the curves and contours of the skin and body, sticking with your patients wherever they go.

Product		Reference Number	Box Contents	Case Contents	Description
Strips		00230 00231 76448000	100 100 300	12 Boxes 12 Boxes 6 Boxes	3/4" x 3" 1" x 3" 1" x 3"
Spots		00301 00303	100 100	12 Boxes 12 Boxes	7/8" Round 1 1/4" Oval
Patches		00330 00340000 78010-00	100 50 50	12 Boxes 12 Boxes 12 Boxes	1 1/2" x 2" 2" x 3" 4" x 2 3/4"
Knuckles		01390	100	12 Boxes	1 1/2" x 3"
Toe Shield (non-sterile)	T	00801	100	12 Boxes	2 1/2" x 1 3/8"
Digits	LH	01306 01307	100 50	12 Boxes 12 Boxes	Sm. Finger Tip Lg. Finger Tip
4-Wing	I	00385	50	12 Boxes	3" x 3"

^{*} Laboratory testing results are intended to illustrate a product's performance under controlled conditions. Actual use results may vary.



Tel. (+1) 704 554 9933; Fax (+1) 704 358 4558 To order toll-free: (+1) 800 552 1157

www.bsnmedical.com

Tel. (+1) 704 554 9933; Fax (+1) 800 835 4325 To order toll-free: (+1) 800 537 1063 www.jobst.com



0100 R13 ©2009 BSN medical Inc. Printed in U.S.A. REV 09/09