

## Innovators for Optimal Lung Care

EasyOne Air is the new, all-in-one revolution in lung function testing. Built for healthcare providers large and small, EasyOne Air puts the power of the most advanced spirometry testing tools in the palm of your hand.

Revolutionary TrueFlow ultrasound technology means unprecedented accuracy for reliable diagnosis and trusted treatment. Equipped with ready-to-use EMR connectivity and a user-friendly interface, the EasyOne Air makes the best in lung care available anywhere.

Learn more about EasyOne Air at nddmed.com 300 Brickstone Square, Suite 604 Andover, MA 01810 | (978) 470-0923 | customerservice@nddmed.com


EasyOne Air

ndd

07.02.2017 15:49



## Simple. Certain. Proven.





The original ultrasonic flow measurement is highly accurate in all flow ranges, independent of gas composition, pressure, temperature and humidity. ndd TrueFlow is a hygienic and resistance-free solution that does not require calibration during its lifetime.



EasyOne Connect software offers one database and one platform for all EasyOne point-of-care-solutions. Data management has never been easier.

## **Specifications**

Test Types FVC, FVL, SVC, MVV and Pre/Post Testing

Languages English, French, German, Spanish

Printing Options Direct to printer or with EasyOne Connect software

EMR connectivity EasyOne Connect

**Dimensions** 3.4x6.1x1.4", 13 oz

Meets ATS/ERS 2005, OSHA and SSD requirements

## Part Numbers

**2500-2A** EasyOne Air spirometry system

**5050-50** EasyOne Flow Tube box of 50

**5050-200** EasyOne Flow Tube box of 200