

Walkabout Four-Wheel Imperial Rollator

The Lumex Walkabout Four-Wheel Imperial (Bariatric) is ideal for larger individuals. Similar to the Walkabout Lite, this rollator offers a comfortable, lightweight and affordable solution to mobility. Designed with an aluminum frame, this lightweight rollator weighs 20 lbs and supports a maximum weight capacity of 450 lb - EVENLY DISTRIBUTED. The Walkabout folds quickly and easily into a compact unit for storage and transport.





Features

- Ergonomic hand grips
- · Easy to operate locking and "quick stop" loop brakes (pull up to stop; push down to lock)
- Padded 20" wide seat for increased seating surface and added comfort
- Removable padded backbar
- Adjustable handle height to accommodate broad range of user heights
- Durable, heavy duty 8" wheel for indoor and outdoor use on even terrain
- Basket included that can be placed underneath seat for carrying of personal items
- Designed to meet the requirements of HCPCS codes: E0149 & E1056
- Includes easy to read, retail-friendly user manual
- Five (5) year limited warranty

And now with added rear reflectors that increase visibility at night!



Straight Backbar

RJ4400K Black 1 ea **RJ4400R** Burgundy 1 ea



Color options

Rollators are not intended to be propelled while seated. Do NOT use this rollator as a wheelchair or transport device.



Specifications

Width Inside Handgrips	. 21.5"
Overall Width	. 28.7"
Overall Depth	26"
measured wheel to	wheel
Seat Width	
Seat Depth	. 13.8"
Seat Height	. 21.8"
Caster Diameter	8"
Minimum Handle Height	34"
Maximum Handle Height	39"
Weight Capacity	450 lb
EVENLY DISTRIE	BUTED
Product Weight	. 20 lb