

LUMEX®

A Graham-Field® Family Brand

Deluxe Forearm Crutches

- Leg and forearm sections can be adjusted independently for a custom fit
- Vinyl-coated arm cuffs can be shaped to the user's arm for extra comfort
- Contoured vinyl handgrip is durable and comfortable
- Telescoping parts have internal bushings and external lock nuts to silence any noise
- Extra-large, non-skid crutch tips provide added traction and stability
- Constructed with high strength aluminum tubing
- Easy push-button height adjustment with locking collet for added security
- Grip-to-floor height adjusts in 1" increments
- Arm cuff to grip height adjusts in 1" increments
- Includes laminated color operating instructions
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: E0110
- Three (3) year limited warranty



6350A	Adult	28" - 37"	1 pr
6350T	Tall	33" - 42"	1 pr
6350Y	Youth	22" - 31"	1 pr

Cuff adjustment from 9" - 12" from handgrip to top of cuff

6350-CUFF	Replacement cuff	1 ea
6350-HG	Replacement handgrips	1 pr
6350-TIP	Replacement tip	1 pr

graham
field™

www.grahamfield.com

© 2017 GF Health Products, Inc.
Atlanta, GA All rights reserved.
GF1700065RevA17