


# Brewer Medical Seating

## Century Series


- Preferred by medical facilities worldwide
- 16" wide, 4" thick poly-foam seat cushion
- Pneumatic height adjustment
- Five-leg, 23" brushed aluminum accent base
- 2" chrome hooded dual wheel casters
- Optional adjustable contoured backrest with black plastic shroud
- Optional adjustable 18" chrome footing



#11001BFR  
in Burgundy

## Value Plus Series


- Affordable without compromising style and reliability
- 15" wide, 4" thick poly-foam seat cushion
- Pneumatic or spinlift height adjustment
- Five-leg, 21" black ABS plastic composite base
- 2" black plastic dual-wheel nylon casters
- Optional adjustable contoured backrest with black plastic shroud



#22500B  
in Dove

## Tradition Series

- Economical exam stools that stand the test of time
- 14" wide, 3.5" thick poly-foam seat cushion
- Sturdy, four-leg, 19.5" chrome base
- 2" rubber wheel, ball bearing casters



#23051  
in Adriatic

# Brewer Medical Seating

## Millennium Series

- 18.5" wide contoured seat cushion
- Adjustable backrest black plastic shroud, excl. 21340 & 21340V
- 2" dual-wheel nylon casters with polyurethane tread; locking casters optional
- 23" five leg powder coated black aluminum base
- Ergonomic design; seamless upholstery standard, excl. 21340

accessories for the millennium series:

#2102265 Set of 5 pressurized locking casters

- B - backrest
- V - seamless/vacuum formed
- A - stationary armrests



#21435BA  
in Azure Blue



#21521B  
in Adriatic



#21340B  
in Agate Blue

### Task

- Pneumatic height adjustment: 18" - 24.75"
- Fully adjustable forward seat tilt
- Adjustable backrest with synchronized traversing movement of seat



#21435BA  
w/ arm rests

#21435B  
w/o arm rests

### Lab

- Pneumatic height adjustment:
  - With Glides: 20.75" - 31"
  - With Casters: 21.75" - 32"
- Standard 18" adjustable foot ring



#21520B w/glides  
#21521B w/casters

### Surgeon

- Hands-free foot pneumatic height adjustment: 19" - 24"
- 18.5" wide contour poly-foam seat cushion on 21340B only; 16" wide on 21340
- Upholstery is stitched on 21340



#21340B  
#21340BV

#21340  
#21340V

The Brewer Company is an ISO 13485:2003 certified company.

Weight capacity for stools - 250 lb.

All Brewer Medical Seating comes with a 5-year warranty and is available in 47 standard colors.



The Brewer Company reserves the right to improve or modify the design of its products consistent with new methods and techniques.

Call your local Brewer sales representative to learn about our complete product line.


The Brewer Company LLC  
N88W13901 Main Street, Menomonee Falls, WI 53051  
1.888.Brewer.1 www.brewercompany.com

# Brewer




## Standard Upholstery Color Selections

For Brewer Tables, Side Chairs, Blood Drawing Chairs and Exam Stools



## California Technical Bulletin 133 Upholstery Color Selections (CAL-133)


For Brewer Tables, Century Series Stools and Value Plus Series Pneumatic Stools



CAL TB 133 is a written legal standard outlining flammability test procedures for upholstered furniture. These vinyl color selections contain a flame blocking barrier that helps to prevent ignition, reflects heat and is flame retardant.

## Premium Plush Upholstery Color Selections

For Brewer Assist Procedure Tables, Access High-Low Exam Tables, Side Chairs and Century Series Stools




All Premium Plush Upholstery Colors are Ultrafabrics ULTRALEATHER™.

**Note:** Tables featuring a plush upholstery top are two-toned, stitched, with Beech and one of the premium plush upholstery colors. Exam stools and side chairs are one solid color, stitched. (-UL)

For Actual Material Swatches call 1•888•Brewer•1

## Additional Standard Upholstery Color Selections

For Brewer Exam Stools Only



For Actual Material Swatches call 1•888•Brewer•1